



Bio-Therapeutic
Skin Care for life inside and out®

bt-titan® MN

advanced system: Model AE

International Edition: Professional System Detail
One Year Limited Warranty

The bt-titan® MN model AE allows for a highly controlled and precise service that is adapted to your clients exact skin needs and comfort. Unlike more invasive procedures, this quick application not only delivers marked results but is suitable for most skin types and tones. With no down time, bt-titan MN applications are a great option for those with a busy lifestyle and schedule.

Unit Specifications:

Power Supply: Rechargeable lithium battery pack or A/C outlet

Power Source: 100-240V; 50/60Hz

Input Power: 5V DC; 1A

Internal Power: 3.7V

Weight: 120g / 4.25oz (with battery)

Dimensions: 155mm long x 21mm round

Sterile and disposable single use cartridges

Quick release cartridge installation

Five variable speeds, Max RPM: 16,000

Length adjustment between 0mm - 1.5mm

Exclusive internal locking mechanism

IP24 water resistant

Power Supply

For charging or external power operating:

AC/DC Adapter model LC-H050100: input 100-240V AC; 50/60Hz;

output: 5V DC; 1A

For internal power:

Rechargeable Lithium battery, model HDY16340; 3.7V DC

Included:

bt-titan MN Model AE

bt-titan MN Operation Manual

bt-titan MN Quick Start Guide

Charging cradle / equipment holder

Disposable cartridges: nano, 12, 36

2 rechargeable lithium battery packs

Antimicrobial and antibacterial silicone carrying case

Protective fabric sleeve

Sanitary plastic sheaths

A/C adapter attachment

Battery charger

A/C adapter

Power cord



Bio-Therapeutic products and equipment are not intended to treat or diagnose any medical conditions.

Items may vary without notice and are not available in all countries.



Bio-Therapeutic
Skin Care for life inside and out®

bt-titan® MN

advanced system: Model MD

International Edition: Professional System Detail
One Year Limited Warranty

The bt-titan® MN model MD allows for a highly controlled and precise service that is adapted to your clients exact skin needs and comfort. Unlike more invasive procedures, this quick application not only delivers marked results but is suitable for most skin types and tones. With no down time, bt-titan MN applications are a great option for those with a busy lifestyle and schedule.

Unit Specifications:

Power Supply: Rechargeable lithium battery pack or A/C outlet

Power Source: 100-240V; 50/60Hz

Input Power: 5V DC; 1A

Internal Power: 3.7V

Dimensions: 155mm long x 21mm round

Weight: 120g / 4.25oz (with battery)

Sterile and disposable single use cartridges

Quick release cartridge installation

Five variable speeds, Max RPM: 16,000

Length adjustment between 0mm - 2.5mm

Exclusive internal locking mechanism

IP24 water resistant

Power Supply

For charging or external power operating:

AC/DC Adapter model LC-H050100: input 100-240V AC; 50/60Hz;

output: 5V DC; 1A

For internal power:

Rechargeable Lithium battery, model HDY16340; 3.7V DC

Included:

bt-titan MN Model MD

bt-titan MN Operation Manual

bt-titan MN Quick Start Guide

Charging cradle / equipment holder

Disposable cartridges: nano, 12, 36

2 rechargeable lithium battery packs

Antimicrobial and antibacterial silicone carrying case

Protective fabric sleeve

Sanitary plastic sheaths

A/C adapter attachment

Battery charger

A/C adapter

Power cord



Bio-Therapeutic products and equipment are not intended to treat or diagnose any medical conditions.

Items may vary without notice and are not available in all countries. The bt-titan MN model MD is not available in the US